


Scan Sound Inc.

1121 S. Military Trail • Suite 118
Deerfield Beach, FL 33442
(954) 425-0199 • FAX (954) 827-2408



Installation Evaluation Form

Please have your M.R.I. installation personnel or Site-Engineer complete the attached evaluation form when placing your order. This data will be used as a guideline in manufacturing the types of cables required, simplifying the Patient Stereo System Installation.

Scan Sound, Inc.

NOTE:


Your valued assistance in completing this form will be greatly appreciated. Answer the questions where appropriate to help Scan Sound determine the best type of installation possible for your facility. If you have any questions, please do not hesitate to call a Scan Sound representative at (954) 425-0199.

Date	Facility	
Department		Phone
Address		Contact Person
City	Order Number	
State	Zip Code	
MRI System Manufacturer		
Person Who Filled Out This Questionnaire		

Please return this form to Scan Sound, Inc.- FAX: (954) 827-2408

Wiring Requirements

The following illustration shows how the wiring is sent into the MRI scan room, via BNC connectors on the RF Penetration panel.



Please refer to the Diagram on the following page when responding to these questions:

1. Does your system currently have two spare BNC connectors that are available on the pass-thru connector panel? YES
 No; We will need two to be shipped with the patient stereo at an additional cost of \$5.95 each
-

2. Is your MRI Magnet Shielded or Unshielded?
 Shielded Unshielded
-

3. What is the field strength of your Magnet
-

4. What length of cable will be required to run from the Stereo Amplifier to the Scan-room Bulkhead Connector Panel?
-

5. What length of cable will be required to run from the Bulkhead Connector panel to the magnet Patient Entry Area?
-

6. What length of cable will be required to run from the AM/FM Receiver to the area where the Microphone is placed at the operator's console?

Choose from the following lengths: 6-feet
 12-feet

7. Will your system require the addition of BNC Connector Pass-Through filters at the Penetration Panel? *Used in MRI Facilities Where RFI Filtering is Recommended by the MRI Manufacturer*
 YES.....Please add a pair to our order and bill us for the added \$299.90 cost
 No

Note: Most systems will not require extra filtering

8. *Who will be designated to install the Stereo System Kits at the facility?
-

*Each site will designate someone to install the system. The installation is well documented and is easily done by anyone familiar with Stereo System hook-up, the ability to run the signal wires via the connector panel and a familiarity with the precautions of working around an energized MRI Magnet.

Scan Sound, Inc.

Please return this form to Scan Sound via Fax at (954) 827-2408

If you prefer to send the form information by email, please send it to: **info@scansound.com**

This form may be filled out online if you prefer.

URL: <http://www.scansound.com/PDF%20Files/EvalForm.pdf>

Cable Measurement Guide

