

**Model J-4002 Stereo
with AUX input from
Computer**

Bose Acoustic Wave
Music System II
AM/FM/CD Player

Speakers Built-In
to Be Heard in
Area Near Console

Music from
Bose Stereo

Aux Input

Computer Used to Stream Music

Mixer / Amplifier
Adjusts Sound for Patient

Microphone
to Speak to
Patient

Shielded Scan Room

RF Penetration
Panel

MRI Magnet

Optional
RF
Filter
Inserts

Optional
Noise
Reduction
Headphone

Headphones
Hung
On Plastic Clips
When Not In Use

Patient
Acoustic
Interface

